

Roche + CDC Partnership (ASLM Webinar)

Develop Capacity and Improve Efficiencies in Laboratory Operations

Asif Ali – Project Manager , Global Access Program (EMEA-LATAM) , Roche Diagnostics

Denise Heaney – Senior Scientific Affairs Manager , Roche Diagnostics (US)

February 11, 2021

A decorative footer featuring a complex geometric pattern of overlapping triangles in various shades of blue, green, pink, and purple. A large, dark blue cross is centered over the pattern.

Purpose and Disclaimer Slide

The purpose of this presentation is to provide information on tools developed in partnership with the CDC ILB to support continued scale up of HIV viral load and early infant diagnosis. These tools provide support training and implementation efforts and can contribute to improving laboratory efficiency and strengthening across the programs.

These tools will be available via the ASLM Resource Center for easier access online providing the ability to access them when it is convenient for your lab, to provide you with updates for the latest innovations to support your HIV programs and to help reduce time spent to build training, implementation or other laboratory activities.


This information is shared for your awareness of their availability, location and content and does not reflect endorsement by ASLM.

Agenda

**Introduction:
Partnerships**

**Implementing cobas[®]
Plasma Separation Card**

**Summary &
Questions**



**Building Capacity on EID DBS to
Enhance Access to testing**


**Streamlining
Verification of Assays**

Agenda

Introduction: Partnerships

Implementing cobas[®]
Plasma Separation Card

Summary &
Questions



Building Capacity on EID DBS to
Enhance Access to testing

Streamlining
Verification of Assays

Public-Private Partnership (PPP) for HIV/AIDS Response

Developing laboratory tools for viral load and EID scale up

Menu Search Media

Roche and the US President's Emergency Plan for AIDS Relief (PEPFAR) partner to strengthen laboratories medicine training and knowledge on the African continent

Basel, 04 December 2012

Today Roche Diagnostics and the US President's Emergency Plan for AIDS Relief (PEPFAR) announced a public-private partnership to improve the capacity of laboratorians and strengthen laboratory systems in African nations highly impacted by HIV/AIDS. Valued at \$12 million over five years, the goal of the partnership is to build institutional capacity of regional and national partners by developing certification programs for laboratory technicians, strengthening pre-service curricula, and institutionalizing quality improvement in diagnostics services.

Shrivastava et al. BMC Health Services Research (2019) 19:46
https://doi.org/10.1186/s12913-018-3744-z

BMC Health Services Research

RESEARCH ARTICLE Open Access

Role of public-private partnerships in achieving UNAIDS HIV treatment targets

Ritu Shrivastava¹, Peter N. Fonjongo^{1*}, Yenew Kebede², Rajendra Bhimaraj³, Shabnam Zavahir³, Christina Mwangi⁴, Renuka Gadde⁵, Heather Alexander¹, Patricia L. Riley¹, Andrea Kim¹ and John N. Nkengasong⁶

Abstract
Background: Despite progress towards achieving UNAIDS 90–90–90 goals, barriers persist in laboratory systems in sub-Saharan Africa (SSA) restricting scale up of early infant diagnosis (EID) and viral load (VL) test monitoring of patients on antiretroviral therapy. If these facilities and system challenges persist, they may undermine recorded gains and appropriate management of patients. The aim of this review is to identify Public Private Partnerships (PPP) in SSA that have resolved systemic barriers within the VL and EID treatment cascade and demonstrated impact in the scale up of VL and EID.

Improved laboratory capacity and accelerated training of lab staff

Quality improvement for laboratory services, including certification courses for molecular diagnostics & QMS



CDC
ASLM

Developing pre-service training curricula

Collaboration with the ASLM to strengthen local capacity and to promote country ownership and sustainability

“PPP contributions to laboratory systems strengthening are a model”¹

Source: <https://www.roche.com/media/releases/med-cor-2012-12-04.htm>. Accessed 5 Sept 2018

1. Shrivastava et al. BMC Health Services Research (2019); <https://doi.org/10.1186/s12913-018-3744-z>

Agenda

Introduction:
Partnerships

**Building Capacity on EID DBS to
Enhance Access to testing**

Implementing cobas[®]
Plasma Separation Card

Streamlining
Verification of Assays

Summary &
Questions



Early Infant Diagnosis Through DBS Collection

Training remains a critical component to expanding access to EID

UNAIDS Joint Statement:

Recently issued by global partners reinforced the commitment to urgently address the challenges to the early diagnosis and treatment of children living with HIV, which led to **estimated death of 95,000 children in 2019¹**.

Coalition of governments and global partners have, in accordance with WHO guidelines, effectively implemented EID across multiple settings, with **training and capacity building of healthcare workers** on DBS Collection across PMTCT sites a key intervention.

Critical success factors to ensuring accurate, timely early diagnosis include the following: i) **obtaining high quality samples**, ii) **careful handling of specimens**; and iii) **maintaining accurate records**.



[1] UNAIDS. [Joint statement calling for urgent country scale-up of access to optimal HIV treatment for infants and children living with HIV](#). 22 December 2020

ASLM Resource Centre

Training material to enable critical successful factors for EID

ASLM Resource Centre Page

The screenshot shows the ASLM Resource Centre website. At the top, there is a navigation bar with links for 'Who We Are', 'What We Do', 'Membership Opportunities', 'Resource Centre', and 'COVID-19'. Below the navigation bar is a banner for the 'Resource Centre' with a date of 'July 31, 2019'. The main content area features a video titled 'Early Infant Diagnosis Through Dried Blood Spot Collection' dated 'June 17, 2019'. Below the title are social media sharing buttons for Facebook, Twitter, and LinkedIn. The video description states: 'Early infant diagnosis (EID) is critical for timely initiation of antiretroviral treatment (ART) in HIV-infected children who are at high risk of mortality. In recognition of the immense benefits of dried blood spot (DBS) as a means of increasing the access to EID, the World Health Organization, African Society for Laboratory Medicine, Centers for Disease Control and Prevention, United States Agency for International Development, Clinton Health Access Initiative, Médecins Sans Frontières, Global Fund, Riders for Health, Clinical Laboratory Standards Institute convened, identified and developed tools that would assist programs in scaling up HIV Diagnosis and uptake of results. This video, produced by the US Centers for Disease Control and Prevention English, French, and Portuguese, describes the procedure for collection of a DBS specimen from an infant below 18 months of age for performing HIV-1 PCR.' Below the text is a video player showing a woman in a red hat holding a baby while a healthcare worker in a white uniform and vest looks on. The video player has 'English' selected, 'DBS EID English' as the title, and 'Watch later' and 'Share' buttons.

- In collaboration with global partners, CDC has led the development of training videos covering **EID through DBS collection**.
- In addition to explaining how to troubleshoot common issues, the videos outline how to:
 - gather DBS supplies for EID*
 - collect DBS samples*
 - tell the difference between valid and invalid DBS; and*
 - ensure DBS cards are properly labelled and logged.*
- A pre- and post-test is available to help track progress
- In English , French and Portuguese, the videos can be accessed via [ASLM's Resources Centre](#) and [YouTube](#).^{1,2}

[1] ASLM. <https://aslm.org/resource/early-infant-diagnosis-through-dried-blood-spot-collection/>. 17 June 2019

[2] ASLM. <https://www.youtube.com/watch?v=2FUlphI5PoU>. 10 May 2018

Agenda


Introduction:
Partnerships

Building Capacity on EID DBS to
Enhance Access to testing

**Implementing cobas[®]
Plasma Separation Card**

Streamlining
Verification of Assays

Summary &
Questions




Limitations to Meeting Virally Suppressed Target - 95% by 2025

Barriers to plasma-quality Viral Load for patients across all settings

Not all people live close to a major health clinic

Infrastructure not in place to allow for optimal storage and transportation of EDTA Plasma Specimens i.e. cold chain



Those living in rural settings have blood taken in small clinics – centrifuge required for sample handling

DBS is an alternative sample type but a more effective solution is required

Use of DBS Can Cause Misclassification

Ideal for EID but not VL

Data from Vietnam shows 6% of patients misclassified using DBS



*Downward
misclassification*

Believe patients failing
treatment **are suppressed**

This could lead to **additional health complications** and the risk of infecting more people which both add to cost

*Upward
misclassification*




Believe suppressed patient
is **failing treatment**

This leads to **additional costs** associated with adherence counselling and the change to more costly second line therapy

cobas[®] Plasma Separation Card

An effective alternative to DBS...



WHO standard of <1000 cp/mL (in plasma)
No other CE-marked or WHO prequalified sample collection solution matches performance



Superior solution compared to DBS
Better performance than DBS with low misclassification rates



Easy handling, transport and storage
Provides an easy-to-use solution to monitor patients in remote areas





Plasma-based sample (gold standard)
*Separates from whole blood
 Stabilizes dried plasma/viral RNA under extreme heat and humidity*



Correlation to plasma as the GOLD standard
Plasma is the preferred sample type for HIV viral load testing


Design of the cobas[®] Plasma Separation Card



Sample Collection Process for the cobas[®] Plasma Separation Card


1) Sample Collection

(140µL whole blood required per spot and one spot required)



2) Sample Transport from the Field to the Laboratory

(18-45°C and up to 85% humidity (>21 days))



cobas[®] Plasma Separation Card Bundle Components

Convenient for clinics and patients, and simplifies ordering



Training and Support for cobas® Plasma Separation Card

Tools to support implementation of cobas® PSC



Roche has developed **training videos, posters and presentations** detailing how to use the **cobas® Plasma Separation Card**, specifically:

- Sample Collection Method
- Preparation for Sample Transportation
- Sample Preparation and Instrument Workflow for Roche systems (including **cobas® 6800/8800 Systems**, **cobas® 4800 System** and **COBAS® AmpliPrep/COBAS® TaqMan® System**)

Training material is available in English , French and Portuguese, and can be accessed via [**ASLM's Resources Centre**](#)

Mozambique Study

Performance - Very low misclassification rates with **cobas**[®] PSC

Misclassification rates of viral load values below 1000 c/ml



0.7%

Capillary specimen




0.3%

Venous specimen



5.9%


DBS specimen



Eliminates the interference of cell-associated viral nucleic acid



No over-quantification



Cell-associated viral nucleic acid present in whole blood



Over-quantification

Milestones in regulatory process met

Evaluation Letter

- **Memorandum issued in 2019 by CDC** following evaluation study in South Africa


WHO Prequalification

- WHO Prequalified for **cobas[®] 6800/8800 Systems** and **COBAS[®] AmpliPrep/COBAS[®] TaqMan[®] System**, application submitted for **cobas[®] 4800 System**
- The **Limit of Detection estimated at 489 cp/mL** (95% CI: 330-1256) for cobas Plasma Separation Card (cobas[®] 6800/8800 Systems)

Verification Panels

- In partnership with **CDC development of verification panels (TBD)**

[1] WHO. [WHO PQ of IVD - cobas HIV-1 Quantitative nucleic acid test for use on the cobas 6800/8800 Systems](#). October 2020



Agenda

Introduction:
Partnerships

Implementing cobas[®]
Plasma Separation Card

Summary &
Questions



Building Capacity on EID DBS to
Enhance Access to testing

**Streamlining
Verification of Assays**

Roche Assay Implementation Tool

Streamlining the verification of new assays for laboratories

Planning assistance for assay verification experiments for PCR based testing for CE-IVD or FDA approved assays*


Data analysis of the performance specifications required for assay verification through a web-based portal

Report generation of analyzed data for assessment by authorized laboratory personal for verification sign-off and compliance requirements

**Tool in current format was developed for Roche PCR based IVD assays*

Roche Assay Implementation Tool

A support tool for laboratories across Africa for verification of IVDs



With the Assay Implementation Tool* labs will be able to:

- **Review** definitions of various verification terms as well as recommended guidelines **on how to perform studies** when verifying the performance of an *In Vitro* Diagnostic test
- **Input data** generated from studies performed as part of a verification process as well as **automatically calculate and produce standard graphs**
- **Document** (via .pdf format) for the laboratory's records the runs were performed as part of a verification study, including the results, data analyses generated, as well as other details associated

* Roche has worked in collaboration with CDC ILB team to modify the tool to meet verification parameters required per CDC guidelines.

Roche Assay Implementation Tool

Accessible, easy to navigate, simple data input interface and report generation

Pathogen & Test Parameters & Definitions

Validation Verification Parameters

Please select parameters for test. For a definition of test, please mouse over the test name.

Limit of Detection
Definition:
 • The lowest concentration of analyte that can be consistently detected (typically, in ≥ 95% of samples tested under routine clinical laboratory conditions) in a defined type of sample.²

How LOD will be performed:
 • A panel member will be chosen from the linearity panel and 20 replicates will be run to generate results for comparison to the manufacturer's reported limit of detection (LOD).⁴

Method Correlation Results

The results from the method correlation are shown using Deming Regression analysis in a graphical format. Along with the graph is also a regression equation using the linearity data so it is possible to assess how close the observed values are to the expected values (Ideally an R² value as close to 1.00 as possible). A line of unity is also included in the graph to show where results that are equivalent would be plotted.

The purpose of the evaluation is to determine if the two methods yield equivalent results within the statistical power of the experiment

If the comparative method is the same as the one used by the manufacturer in the statement of claims, it is possible to compare statistically the experimental results to the manufacturer's claim to verify acceptable performance.¹

A minimum of 40 patient samples with

Method Correlation Quantitative Assessment - cobas® HIV-1 for use on the cobas® 6800 System

Pathogen & Test Parameters & Definitions

Pathogen & Test

- HBV Hepatitis B Virus
- HCV Hepatitis C Virus
- HIV-1 Human Immunodeficiency Virus
- HPV Human Papillomavirus

- COBAS® AmpliPrep/COBAS® TaqMan® HCV Test, v2.0
- cobas® HCV for use on the cobas® 4800 System
- cobas® HCV for use on the cobas® 6800 System
- cobas®

Homepage

Please note: This tool is to be used as a guide and recommendations do not serve as definitive ways in which a laboratory should perform validation studies or use the analyses performed. Local guidelines and recommendations, where available, should supersede anything that is included in this tool.

This tool is provided by (the CDC and/or ASLM) as a reference for the laboratory's use in conducting assay evaluation and verification. The laboratory's use of data collected by this tool should not be considered a substitute for the laboratory's own independent review and verification of the assay's performance, and does not ensure compliance with applicable laws, regulations, accreditation requirements or international standards; such compliance remains solely the responsibility of the laboratory. CDC, ASLM and Roche Molecular Systems, Inc. expressly disclaim any liability resulting from the laboratory's use or approval of this tool or any data generated by this tool.

Welcome to the Validation Verification Tool Website

This tool is designed to help laboratories with implementing new tests making it easier to design, setup, and analyze data generated as part of a new test validation verification procedure.

When a laboratory decides to implement or evaluate a new molecular diagnostic test, it is important to be able to perform several validation verification studies. The decision on how to approach this process resides with the laboratory director and is based upon local regulatory and legal guidelines. This tool is intended to help supplement the validation verification of a new test.

Getting Started

The Validation Verification tool is intended to help laboratories define, setup, and test various validation parameters as part of the process for verifying the performance of an *In Vitro* Diagnostic test.

By using this tool, laboratory users will be able to:

- Review definitions of various validation terms as well as recommended guidelines on how to perform validation studies when verifying the performance of an *In Vitro* Diagnostic test.
- Be able to input data generated from studies performed as part of a validation verification process as well as automatically calculate and produce standard graphs.

Jump to

- Select Pathogen & Test for Validation Verification
- Edit/View Archived Tests
- Quantitative Tests: Parameters & Definitions
- Qualitative Tests: Parameters & Definitions

Precision Results

The laboratory is responsible for establishing the precision of each system by assessing day-to-day, run-to-run and within-run variation, as well as operator variance.

Accomplished by repeat testing of known patient samples over time, testing QC material over time and/or repeat testing of calibration materials over time.¹

Two levels (highest and mid-range titer) of diluted control materials plus negative control can be run in triplicate one run per day for 5 days.

Total precision is assessed as the percentage of coefficient of variability; the standard deviation expressed as a percentage of the mean value of the replicate measurements.

The precision can be established within the run, run to run, and as a total precision.^{1,2}

Manufacturer's precision data can be found in the package insert under the section: Non-Clinical Performance Evaluation

Precision Calculations - cobas® HIV-1 for use on the cobas® 6800 System

Precision data comparison

No	Sample	Expected Result, copies/mL	Log ₁₀ Expected Result	Observed Result, copies/mL	Log ₁₀ Obs Result	Diff. Log ₁₀
1	D1 DIL OL3 R1	1,000	3	1,760	3.25	0.25
2	D1 DIL OL3 R2	1,000	3	1,460	3.16	0.16
3	D1 DIL OL3 R3	1,000	3	979	2.99	-0.01
4	D4 DIL OL3 R1	1,000	3	1,350	3.13	0.13
5	D4 DIL OL3 R2	1,000	3	2,150	3.33	0.33
6	D4 DIL OL3 R3	1,000	3	2,160	3.33	0.33
7	D1 DIL OL6 R1	1,000,000	6	2,210,000	6.34	0.34
8	D1 DIL OL6 R2	1,000,000	6	1,480,000	6.17	0.17
9	D1 DIL OL6 R3	1,000,000	6	1,540,000	6.19	0.19
10	D4 DIL OL6 R1	1,000,000	6	1,780,000	6.25	0.25
11	D4 DIL OL6 R2	1,000,000	6	1,710,000	6.23	0.23
12	D4 DIL OL6 R3	1,000,000	6	1,810,000	6.26	0.26
13	D2 DIL OL3 R1	1,000	3	2,250	3.35	0.35
14	D2 DIL OL3 R2	1,000	3	1,690	3.23	0.23
15	D2 DIL OL3 R3	1,000	3	1,850	3.27	0.27
16	D5 DIL OL3 R1	1,000	3	2,090	3.32	0.32
17	D5 DIL OL3 R2	1,000	3	2,120	3.33	0.33

Notes:

- Training video on using the tool (which will provide a technical overview) will be issued in parallel to the tool's release

Agenda


Introduction:
Partnerships

Building Capacity on EID DBS to
Enhance Access to testing

Implementing cobas[®]
Plasma Separation Card

Streamlining
Verification of Assays

Summary &
Questions



Summary & Questions

Tools to support capacity development and improve efficiencies in laboratory operations

Capacity building training materials available on ASLM Resource Center

- **DBS collection for EID testing**
- **cobas® Plasma Separation Card** (dried plasma spots), sample collection, transport and laboratory process and country implementation

Introducing a tool currently in development to help optimize laboratory operations (*timing for release to be determined*)

- Online, intuitive **Roche Assay Implementation Tool** to streamline planning and implementation of the verification process

Thank you to the contributors

CDC



Clement Zeh

Role: Team Lead HIV Viral Load & EID| International Lab Branch

Email: cbz2@cdc.gov



Monte Martin

Role: Lead, Laboratory Systems Strengthening Team

Email: mdmartin@cdc.gov

Roche



Asif Ali

Role: Global Access Program Project Manager, SSA

Email: asif.ali.aa4@roche.com



Denise Heaney

Role: Senior Scientific Affairs Manager, Medical & Scientific Affairs

Email: denise.heaney@roche.com



Joni Zurawinski

Role: Director, Access & Partnerships, LMICs EMEA LATAM

Email: joni.zurawinski@roche.com

References

- | | | | |
|----------|--|---|------------------------------|
| <i>a</i> | Shrivastava R et al. (2019) | Role of public-private partnerships in achieving UNAIDS HIV treatment targets | BMC Health Services Research |
| <i>b</i> | Roche and the US President's Emergency Plan for AIDS Relief (PEPFAR) partner to strengthen laboratories medicine training and knowledge on the African continent | From https://www.roche.com/media/releases/med-cor-2012-12-04.htm Accessed 5 Sept 2018 | Roche |
| <i>c</i> | Joint statement calling for urgent country scale-up of access to optimal HIV treatment for infants and children living with HIV | From https://www.who.int/news/item/22-12-2020-joint-statement-calling-for-urgent-country-scale-up-of-access-to-optimal-hiv-treatment-for-infants-and-children-living-with-hiv
Accessed 22 December 2020 | WHO |
| <i>d</i> | Vubil A et al. (2020) | Accurate HIV viral load measurement in primary health care settings using the cobas® plasma separation card. | PLoS ONE 15(5): e0232122. |
| <i>e</i> | WHO PQ of IVD - cobas HIV-1 Quantitative nucleic acid test for use on the cobas 6800/8800 Systems | https://www.who.int/diagnostics_laboratory/evaluations/pq-list/hiv-vrl/201026_amended_pqpr_0365_118_00_cobas_ .
Accessed October 2020 | WHO |

Doing now what patients need next